



Disclosures

- STERIS Corporation is providing the speakers and continuing education credits for this activity. Presenters are employees of STERIS Corporation and receive no direct compensation other than their normal salaries for participation in this activity.
- Commercial products referred to or seen during this presentation do not constitute a commercial support by the speakers.

Continuing Education

- STERIS Corporation is an approved provider of continuing nursing education by **CBRN** – provider # CEP 11681 and an approved Administrator Education Unit (AEU) and Infection Prevention Control (IPCH) provider by **BASC** – provider # 1417.
- This program is approved for:
 - **0** hour(s) of GI Specific content credit by **ABCGN** (American Board of Certification for Gastroenterology Nurses),
 - **1** AEU(s) & **1** IPCH(s) by **BASC** (Board of Ambulatory Surgery Certification), and
 - **1** contact hour(s) of continuing education credit
 - **CBRN** (California Board of Registered Nursing);
 - **CBSPD** (Certified Board for Sterile Processing and Distribution); and
 - **IAHCSMM** (International Association of Healthcare Center Service Material Management).
- Participants must be present for the entire presentation/seminar to receive continuing education credit; partial credit will not be given.

Partnerships



Through a partnership with CCI[®], it also meets CNOR[®] and CSSM[®] recertification requirements for perioperative nurses.



Learning Objectives

- Identify sterilization failures
- Identify preventive actions to prevent sterilization failures
- Identify how sterilization failures fit into the departmental Quality Management System (QMS) program

Sterilization Modalities

- High temperature
 - Steam
- Low temperature
 - Vaporized hydrogen peroxide
 - Vaporized hydrogen peroxide with gas plasma
 - Vaporized hydrogen peroxide with ozone
 - Ethylene oxide
 - Ozone





What happened?

- Investigation initiated
- Equipment checked – all good here
- Assembly checked – nothing out of the ordinary
- Hmm, can't figure it out
- I hope it doesn't happen again

Reactive Outcomes

- Rushed actions
- Time consuming
- Root cause not identified & fixed
- Customer dissatisfaction
- Negative work environment
- Staff dissatisfaction
- Potential negative patient outcomes



QMS in a Healthcare Organization


ANSI/AAMI ST90:2017

"Minimum Requirements to effectively, efficiently, and consistently process medical devices to prevent adverse patient events and non-manufacturer-related device failure."

"Minimum Requirements to effectively, efficiently, and consistently process medical devices to prevent adverse patient events and non-manufacturer-related device failure."

QMS Language

- Instructions for use
- Non-conformity
- Corrective action
- Preventive action
- Risk assessment

A red book titled "DICTIONARY OF THE ENGLISH LANGUAGE" is shown. The book has a red cover and a spiral binding on the right side. The title is printed in white capital letters on the front cover. The book is open, showing the pages.

-
- A red hardcover dictionary titled "DICTIONARY OF THE ENGLISH LANGUAGE". The book is shown from a three-quarter perspective, highlighting its red cover and the spiral binding on the right edge. The pages are visible, showing some text and a spiral binding structure.

Instructions for Use (IFU)

- Written recommendations provided by the manufacturer of a device that provide instructions for operation & safe & effective processing



DISPOSABLE AND REUSABLE INSTRUMENTS THAT MUST BE USED FOR DISINFECTING SURFACES (This IFU applies to all disinfectant products registered under EPA registration number 7689-22-1010.)

Disinfectant	Preparation	Application	Precautions for Handling	Additional Labeling
Disinfectant	1. The product is supplied as a ready-to-use solution. 2. It contains 1% available chlorine.	1. Apply the product to the surface to be disinfected. 2. Allow the product to remain on the surface for at least 10 minutes. 3. Rinse the surface with water after application.	1. Avoid contact with eyes, nose, mouth, or skin. 2. If contact occurs, flush with water immediately. 3. Do not ingest the product.	
Disinfectant	1. The product is supplied as a concentrate. 2. It contains 1% available chlorine.	1. Dilute the product with water according to the label directions. 2. Apply the diluted product to the surface to be disinfected. 3. Allow the product to remain on the surface for at least 10 minutes. 4. Rinse the surface with water after application.	1. Avoid contact with eyes, nose, mouth, or skin. 2. If contact occurs, flush with water immediately. 3. Do not ingest the product.	
Precautions for Handling	1. Do not mix the product with other chemicals, especially bleach, ammonia, acids, alcohols, or petroleum solvents. 2. Do not use the product on surfaces that are already wet or greasy. 3. Do not use the product on surfaces that are damaged or peeling. 4. Do not use the product on surfaces that are made of metal, stone, or wood. 5. Do not use the product on surfaces that are painted or varnished. 6. Do not use the product on surfaces that are made of plastic or rubber. 7. Do not use the product on surfaces that are made of fabric or paper. 8. Do not use the product on surfaces that are made of leather or fur. 9. Do not use the product on surfaces that are made of glass or crystal. 10. Do not use the product on surfaces that are made of brass or copper.			
Additional Labeling	1. Read the entire label before using the product. 2. Follow the directions on the label carefully. 3. Do not use the product if you are pregnant, nursing, or have a medical condition. 4. Do not use the product if you are allergic to chlorine or bleach. 5. Do not use the product if you have asthma or other respiratory conditions. 6. Do not use the product if you have sensitive skin or allergies. 7. Do not use the product if you have open wounds or cuts. 8. Do not use the product if you have a fever or flu-like symptoms. 9. Do not use the product if you are taking medication. 10. Do not use the product if you are drinking alcohol.			

- [illegible]

Nonconformity

- The result of nonfulfillment of a specified requirement




-
- accept only if negative circle turns orange/yellow
- Chemical Indicator**
- CE
- ISO 9001
- SERIAL 036370
- Made in Italy
- 2020-11-01




Corrective Action


- Action taken to eliminate the cause of a detected nonconformity or other undesirable situation



Dirty device



Correction




Corrective action

-
- A close-up photograph of a metal fastener, likely a bolt or screw, showing significant corrosion. The threads and the surface of the fastener are heavily encrusted with a brown, flaky material, indicating severe rust or degradation. The background is a light blue surface.



```
graph LR; subgraph Causes [Possible Causes]; M[Materials]; E[Environment]; P[Process]; U[Utilities]; PR[Personnel]; EQ[Equipment]; end; M --> Main; E --> Main; P --> Main; U --> Main; PR --> Main; EQ --> Main; Main --> Effect{{Sterilization Failure}}
```



Preventive Action

- Action to eliminate the cause of a potential nonconformity or other undesirable potential situation




Clogged spray arm May give dirty devices PM of spray arms

Preventive Action Outcomes

- Standardized actions
- Scheduled actions
- Preventive measures



Corrective Vs. Preventative Actions



Prevention = Positive Outcomes



Where Do We Begin?



Risk

The possibility of a non-conformity occurring that may lead to harm

Risk Management

Identification and mitigation of events that may lead to harm



Develop the Team

- Sterile processing staff
- Peri-operative staff
- Infection prevention staff
- Risk management/quality management
- Facilities/biomed
- Vendor/manufacturer (ad hoc)

The Meeting

- Start-up meeting
- Team goal
- Meeting schedule
- Expectations for;
 - Attendance
 - Participation
 - Reporting



Identify Non-Conformity/ Harm

- Process and procedures
- Staff activities
- Department feedback, reports and inspections
- Customer complaints
- Recalls and incident reports
- Quality/infection prevention

Root Cause Analysis

- Event or situation that creates the nonconformity
- Can have multiple root causes









Root Cause Analysis Methods




- 5 Whys
- Cause & Effect Diagram/Ishikawa Diagram (Fishbone)
- Failure Modes & Effects Analysis (FMEA)
- Hazard Analysis & Critical Control Point (HACCP)
- Fault Tree Analysis
- Additional information: <http://asq.org>

Process Mapping



Process Flow Chart










Failure Mode	Failure Cause	Failure Effect	
		Local Effect	End Effect
Residual soil on device following cleaning	Improper loading of washer - disinfectant prevents soil removal	Device must be recleaned	Added labor and cost Delayed procedure
		Residual soil prevents sterilization of the device	Cross contamination resulting in infection Cross contamination resulting in infection & death
	Sprayer arm blocked	Device must be recleaned	Added labor and cost Delayed procedure
		Residual soil prevents sterilization of the device.	Cross contamination resulting in infection Cross contamination resulting in death





Failure Mode	Failure Cause	Failure Effect		Current Risk Control	Current Risk Level		
		Local Effect	End Effect		S	P	
Residual soil on device following cleaning	Improper loading of washer-disinfector prevents soil removal	Device must be re-cleaned	Added labor and cost Delayed procedure		N	HL	3
					M	HL	2
					CR	NL	3
		Residual soil prevents sterilization of the device.	Cross contamination resulting in infection • Magnified inspection • Prophylactic antibiotics given		CA	NL	4
	Sprocket arm blocked	Device must be re-cleaned	Added labor and cost Delayed procedure		N	HL	1
					M	NL	2
		Residual soil prevents sterilization of the device.	Cross contamination resulting in infection • Magnified inspection		CR	NL	3
		Residual soil prevents sterilization of the device.	Cross contamination resulting in infection & death • Prophylactic antibiotics given		CA	NL	4
	Did not follow all steps in endoscope IFU	Device must be re-cleaned	Added labor and cost Delayed procedure		N	NL	1
		Residual soil prevents sterilization of the device.	Cross contamination resulting in infection • Visual inspection		CR	HL	3
		Residual soil prevents sterilization of the device.	Cross contamination resulting in infection & death • Visual inspection		CA	NL	4



Order of Mitigation

1. Remove the risk

- Reusable to disposable
- Ethylene oxide to vaporized hydrogen peroxide


2. Minimize the probability of occurrence

- Preventative maintenance
- QC tests and inspections


3. Warn of residual risk

- Procedures
- Placards

Examples of Risk Mitigation



Examples of Risk Mitigation



Failure Mode	Failure Cause	Failure Effect		Current Risk Control	Current Risk Level			Proposed Mitigation	New Risk Level		
		Local Effect	End Effect		S	P	Risk		S	P	Risk
Residual soil on device following cleaning	Improper loading of washer - disinfectant prevents soil removal	Device must be reloaded	Added labor and cost Delayed procedure	• Magnified inspection	M	HL	3	• Magnified inspection • Protein detection test • Bore-scope inspection • Prophylactic antibiotics given	M	HL	3
		Residual soil prevents sterilization of the device.	Cross contamination resulting in infection	• Magnified inspection	CR	NL	3		CR	1	2
			Cross contamination resulting in infection & death	• Prophylactic antibiotics given	CA	NL	4		CA	1	3
		Device must be reloaded	Added labor and cost Delayed procedure	• Magnified inspection	N	L	3		N	L	3
	Sprayer arm blocked	Residual soil prevents sterilization of the device.	Cross contamination resulting in infection	• Magnified inspection	CR	NL	3	• Protein detection test • Bore-scope inspection • Prophylactic antibiotics given	CR	1	3
		Residual soil prevents sterilization of the device.	Cross contamination resulting in infection & death	• Prophylactic antibiotics given	CA	NL	4		CA	1	3
		Device must be reloaded	Added labor and cost Delayed procedure	• Magnified inspection	N	NL	1		N	HL	3
		Residual soil prevents sterilization of the device.	Cross contamination resulting in infection	• Visual inspection	CR	HL	3		CR	NL	3
	Did not follow all steps in endoscope IFU	Residual soil prevents sterilization of the device.	Cross contamination resulting in infection & death	• Visual inspection	CA	NL	4	• Visual inspection • Protein detection test • Bore-scope inspection	CA	NL	4
		Device must be reloaded	Added labor and cost Delayed procedure	• Visual inspection	M	NL	2		M	L	2
		Residual soil prevents sterilization of the device.	Cross contamination resulting in infection	• Visual inspection	CR	HL	3		CR	NL	3
		Residual soil prevents sterilization of the device.	Cross contamination resulting in infection & death	• Visual inspection	CA	NL	4		CA	NL	4

Closing the Gap

- Identify the gaps
- Prioritize
- Plan for implementation

Monitor

- Effective mitigation
- Effective corrective action
- Confirms risk assessment assumptions
- Observe for new risks or changes in non-conformity frequencies

The purpose of risk management is to control risk




The purpose of monitoring is to measure the effectiveness of controls



Monitor

- Schedule dependent on risk assessment team's recommendations
- Scheduled reviews of
 - The risk management program
 - Risk assessments
 - Individual risks
 - Corrective actions
- Update risk documents, as needed
- Report to designated departments

Chaos Becomes Order with QMS



Stop Reactivism. Move to Preventionism

- Read your facility's risk management policies and procedures
- Identify your team
- Kick-off your sterilization failure risk assessment

Other Resources



References

- Association for the Advancement of Medical Instrumentation (AAMI), Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities, ANSI/AAMI ST79:2017.
- Association for the Advancement of Medical Instrumentation (AAMI), Processing of Healthcare Products – Quality Management Systems for Processing in Health Care Facilities, ANSI/AAMI ST90:2017.
- Occupational Safety and Health Administration (OSHA). Medical & First Aid Standard. Code of Federal Regulations, Title 29, Part 1910.133
- Occupational Safety and Health Administration (OSHA). Eye and Face Protection Standard. Code of Federal Regulation, Title 29, Part 1910.151
- AORN 2018 Guideline for Sterilization
- <http://www.haccpalliance.org/sub/haccpmodels/guidebook.pdf>

Questions



EDUCATION

STERIS
UNIVERSITY



university.steris.com

Join thousands of your colleagues who earn over
80,000 CE hours each year with STERIS University!

One Integrated Approach to Healthcare.

STERIS